



INNOVATION FOR **SAFE** GYN SURGERY



LapBox[®] POWER


ADVANCED DUAL-WALL CONTAINMENT SYSTEM FOR POWER MORCELLATION



Protecting Patients from Cancer Spillage with Dual-Layer Design

- + Accommodates FDA safety warnings
- + Fast tissue encapsulation into the chamber
- + Ensure **top-quality patient care**
- + **Peace of mind** in every procedure

LapBox[®] is a **next-generation containment system** designed to enable **safe, more advanced laparoscopic GYN surgeries**. With **FDA clearance for power morcellation**, its **dual-wall, airbag-like design** delivers **protection, precision, and confidence** – without compromising workflow.




Pressure

Bag in Bag



For Power Morcellation

LapBox POWER



Required Tools:

5 mm × 30° Laparoscope
10 mm Laparoscopic Claw Grasper

Works only in the umbilical incision

Required training before use

FDA-Cleared for Power Morcellation



Dual-Wall, Spill-Proof Chamber



Inflatable Air-Lock for Organ Protection



Blade-Safe Design fixed-size chamber



Intuitive, easy-to-control delivery system

LEAD THE SHIFT TOWARD SAFE LAPAROSCOPIC TECHNIQUES

Distributor:

LiNA Medical USA
(855) 546-2633 | orders@linamed.com

Part # Ark-P01

Unit of Measure: 5 kits per box



[View Surgical Demo](#)

"LapBox provides the confidence we need to perform the safest and most efficient tissue removal possible"

Amy Garcia, MD, GYN Surgeon,
University of NM Health Sciences Center